# Tracking Bracket

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
**INTENDED USE**

- The device is intended to provide physicians with tools for electromagnetic tracking of instruments with respect to image data.

**CAUTION**

- *Federal (United States) law restricts this device to sale by or on the order of a physician.*

**WARNING**

- *Transducer covers containing natural rubber latex may cause allergic reactions. In the USA, refer to FDA Medical Alert MDA91-I.*
- *Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's manual.*
- *Disposable components are packaged sterile and are single-use only. Do not use if integrity of packaging is violated or if expiration date has passed.*
- *Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.*

**I. Attaching Bracket to Transducer****WARNING**

- *If bracket does not fit properly, do not use.*
- *Ensure bracket is placed on transducer prior to covering.*

**PRIOR TO USE**

- *Refer to your system's manual for use of general purpose sensor.*

1. Attach bracket to transducer by aligning locating features. Ensure bracket is firmly attached.

**II. Attaching General Purpose Sensor to Bracket****WARNING**

- *Ensure device is securely attached, failure to do so could result in inaccuracy.*

1. Place general purpose sensor onto sensor mount area.

**III. Covering Transducer, General Purpose Sensor and Bracket****WARNING**

- *Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.*
- *For illustration purposes only, components may be shown without a cover. Always place a cover over transducer, general purpose sensor and bracket to protect patients and users from cross-contamination.*

1. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
2. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
3. Inspect cover to ensure there are no holes or tears.

**IV. Removing General Purpose Sensor from Bracket**

1. Remove general purpose sensor from sensor mount area.

**V. Disposal****WARNING**

- *Dispose of single-use components as infectious waste.*
- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*
- *See your system's manual for reprocessing transducer and general purpose sensor between use.*

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COPYRIGHT © 2014 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
MPS GmbH